Medical Coverage Policy | Noninvasive Techniques for the Evaluation and Monitoring of Patients with Chronic Liver Disease



EFFECTIVE DATE: 01 | 01 | 2019 **POLICY LAST UPDATED:** 12 | 04 | 2018

OVERVIEW

Noninvasive techniques to monitor liver fibrosis are being investigated as alternatives to liver biopsy in patients with chronic liver disease. This policy documents coverage for transient elastography imaging for patients with chronic liver disease.

MEDICAL CRITERIA

Not applicable

PRIOR AUTHORIZATION

Not applicable

POLICY STATEMENT

BlueCHiP for Medicare

Transient elastography (FibroScan) imaging is covered for the evaluation of patients with chronic liver disease and it is not covered for monitoring of patients with chronic liver disease, as the evidence is insufficient to determine the effects of the technology on health outcomes.

Commercial Products

Transient elastography (FibroScan) imaging is medically necessary for the evaluation of patients with chronic liver disease and it is not medically necessary for monitoring of patients with chronic liver disease, as the evidence is insufficient to determine the effects of the technology on health outcomes.

COVERAGE

Benefits may vary between groups and contracts. Please refer to the appropriate Benefit Booklet, Evidence of Coverage, or Subscriber Agreement for applicable not medically necessary/not covered benefits/coverage.

BACKGROUND

BIOPSY FOR CHRONIC LIVER DISEASE

The diagnosis of non-neoplastic liver disease is often made from needle biopsy samples. In addition to establishing a disease etiology, liver biopsy can determine the degree of inflammation present and can stage the degree of fibrosis. The degree of inflammation and fibrosis may be assessed by different scoring schemes. Most of these scoring schemes grade inflammation from 0 to 4 (0 = no or minimal inflammation, 4 = severe) and fibrosis from 0 to 4 (0 = no fibrosis, 4 = cirrhosis). There are several limitations to liver biopsy, including its invasive nature, small tissue sample size, and subjective grading system. Regarding small tissue sample size, liver fibrosis can be patchy and thus missed on a biopsy sample, which includes only 0.002% of the liver tissue. A noninvasive alternative to liver biopsy would be particularly helpful, both to initially assess patients and then as a monitoring tool to assess response to therapy.

HEPATITIS C VIRUS

Infection with hepatitis C virus (HCV) can lead to permanent liver damage. Before noninvasive tests were available, liver biopsy is typically recommended before the initiation of antiviral therapy. Repeat biopsies may be performed to monitor fibrosis progression. Liver biopsies are analyzed according to a histologic scoring

system; the most commonly used one for hepatitis C is the Metavir scoring system, which scores the presence and degree of inflammatory activity and fibrosis. The fibrosis is graded from F0 to F4, with a Metavir score of F0 signifying no fibrosis and F4 signifying cirrhosis (which is defined as the presence throughout the liver of fibrous septa that subdivide the liver parenchyma into nodules and represents the final and irreversible form of disease). The stage of fibrosis is the most important single predictor of morbidity and mortality in patients with hepatitis C. Biopsies for hepatitis C are also evaluated according to the degree of inflammation present, referred to as the grade or activity level. For example, the Metavir system includes scores for necroinflammatory activity ranging from A0 to A3 (A0 = no activity, A1 = minimal activity, A2 = moderate activity, A3 = severe activity).

HEPATITIS B VIRUS

Most people who become infected with hepatitis B virus recover fully, but a small portion will develop chronic hepatitis B virus, which can lead to permanent liver damage. As with HCV, identification of liver fibrosis is needed to determine timing and management of treatment, and liver biopsy is the criterion standard for staging fibrosis. The grading of fibrosis in hepatitis B virus also uses the Metavir system.

ALCOHOLIC LIVER DISEASE

Alcoholic liver disease (ALD) is the leading cause of liver disease in most Western countries. Histologic features of ALD usually include steatosis, alcoholic steatohepatitis (ASH), hepatocyte necrosis, Mallory bodies (tangled proteins seen in degenerating hepatocytes), a large polymorphonuclear inflammatory infiltrate, and, with continued alcohol abuse, fibrosis and possibly cirrhosis. The grading of fibrosis is similar to the scoring system used in hepatitis C. The commonly used Laënnec scoring system uses grades 0 to 4, with 4 being cirrhosis.

NONALCOHOLIC FATTY LIVER DISEASE

Nonalcoholic fatty liver disease (NAFLD) is defined as a condition that pathologically resembles ALD but occurs in patients who are not heavy users of alcohol. Moreover, NAFLD may be associated with a variety of conditions, including obesity, diabetes, and dyslipidemia. The characteristic feature of NAFLD is steatosis. At the benign end of the disease spectrum, there is usually no appreciable inflammation, hepatocyte death, or fibrosis. In contrast, nonalcoholic steatohepatitis (NASH), which shows overlapping histologic features with ALD, is an intermediate form of liver damage, and liver biopsy may show steatosis, Mallory bodies, focal inflammation, and degenerating hepatocytes. NASH can progress to fibrosis and cirrhosis. A variety of histologic scoring systems have been used to evaluate NAFLD. The NAFLD Activity Score system for NASH includes scores for steatosis (0-3), lobular inflammation (0-3), and ballooning (0-2). Cases with scores of 5 or greater are considered NASH, while cases with scores of 3 and 4 are considered borderline (probable or possible) NASH. The grading of fibrosis is similar to the scoring system used in hepatitis C. The commonly used Laënnec scoring system uses grades 0 to 4, with 4 being cirrhosis.

Transient elastography (FibroScan) uses a mechanical vibrator to produce mild amplitude and low-frequency (50 Hz) waves, inducing an elastic shear wave that propagates throughout the liver. Ultrasound (US) tracks the wave, measuring its speed in kilopascals, which correlates with liver stiffness. Increases in liver fibrosis also increase liver stiffness and resistance of liver blood flow. Transient elastography does not perform as well in patients with ascites, higher body mass index, or narrow intercostal margins. Although FibroScan may be used to measure fibrosis (unlike liver biopsy), it does not provide information on necroinflammatory activity and steatosis, nor is it accurate during acute hepatitis or hepatitis exacerbations.

Transient elastography (FibroScan) is the most widely validated of the noninvasive methods. FibroScan has been studied in populations with viral hepatitis, (NAFLD), and ALD. FibroScan validation studies have suggested that it can provide good detection of significant fibrosis and good-to-excellent detection of cirrhosis compared to liver biopsy for HCV and HBV. There are limited data on NAFLD and ALD. There are no established or validated cutoffs, and the quality of the validation studies was generally not high.

Failures of the test are not uncommon, particularly for those with high BMI; however, failures were frequently missed in analyses of the validation studies. Newer more sensitive probes may lessen this limitation. There is no direct evidence that FibroScan improves health outcomes. However, FibroScan has been allowed as an alternative to biopsy to diagnose fibrosis or cirrhosis to establish trial eligibility in several RCTs that established efficacy of HCV treatments.

For individuals who have chronic liver disease who receive transient elastography, the evidence includes many systematic reviews of more than 50 observational studies (>10,000 patients). Relevant outcomes are test accuracy and validity, morbid events, and treatment-related morbidity. Transient elastography (FibroScan) has been studied in populations with viral hepatitis, nonalcoholic fatty liver disease, and alcoholic liver disease. There are varying cutoffs for positivity. Failures of the test are not uncommon, particularly for those with high body mass index, but these failures often went undetected in analyses of the validation studies. Given these limitations and the imperfect reference standard, it can be difficult to interpret performance characteristics. However, for the purposes of deciding whether a patient has severe fibrosis or cirrhosis, the FibroScan results provide data sufficiently useful to determine therapy. In fact, FibroScan has been used as an alternative to biopsy to establish eligibility regarding the presence of fibrosis or cirrhosis in the participants of several RCTs. These RCTs showed the efficacy of hepatitis C virus treatments, which in turn demonstrated that the test can identify patients who would benefit from therapy. The evidence is sufficient to determine that the technology results in a meaningful improvement in the net health outcome.

CODING

BlueCHiP for Medicare and Commercial Products

The following codes are covered with one of the diagnosis codes indicated below:

76981 Ultrasound, Elastography; Parenchyma (eg Organ) (New code effective 1/1/2019)

91200 Liver elastography, mechanically induced shear wave (eg, vibration), without imaging, with interpretation and report

Diagnosis Codes:

K70.0 through K77

RELATED POLICIES

None

PUBLISHED

Provider Update, January 2019

REFERENCES:

- 1. Crossan C, Tsochatzis EA, Longworth L, et al. Cost-effectiveness of non-invasive methods for assessment and monitoring of liver fibrosis and cirrhosis in patients with chronic liver disease: systematic review and economic evaluation. *Health Technol Assess.* Jan 2015;19(9):1-409, v-vi. PMID 25633908
- 2. Lichtinghagen R, Bahr MJ. Noninvasive diagnosis of fibrosis in chronic liver disease. *Expert Rev Mol Diagn*. Sep 2004;4(5):715-726. PMID 15347264
- 3. Mohamadnejad M, Montazeri G, Fazlollahi A, et al. Noninvasive markers of liver fibrosis and inflammation in chronic hepatitis B-virus related liver disease. *Am J Gastroenterol.* Nov 2006;101(11):2537-2545. PMID 17029616
- 4. Zarski JP, Sturm N, Guechot J, et al. Comparison of nine blood tests and transient elastography for liver fibrosis in chronic hepatitis C: the ANRS HCEP-23 study. *J Hepatol.* Jan 2012;56(1):55-62. PMID 21781944 5. Fraquelli M, Rigamonti C, Casazza G, et al. Reproducibility of transient elastography in the evaluation of liver fibrosis in patients with chronic liver disease. *Gut.* Jul 2007;56(7):968-973. PMID 17255218 Chon YE, Choi EH, Song KJ, et al. Performance of transient elastography for the staging of liver fibrosis in patients with chronic hepatitis B: a meta-analysis. *PLoS One.* Oct 2012;7(9):e44930. PMID 23049764

- 6. Friedrich-Rust M, Ong MF, Martens S, et al. Performance of transient elastography for the staging of liver fibrosis: a meta-analysis. *Gastroenterology*. Apr 2008;134(4):960-974. PMID 18395077
- 7. Kwok R, Tse YK, Wong GL, et al. Systematic review with meta-analysis: non-invasive assessment of non-alcoholic fatty liver disease--the role of transient elastography and plasma cytokeratin-18 fragments. *Aliment Pharmacol Ther.* Feb 2014;39(3):254-269. PMID 24308774
- 8. Shaheen AA, Wan AF, Myers RP. FibroTest and FibroScan for the prediction of hepatitis C-related fibrosis: a systematic review of diagnostic test accuracy. *Am J Gastroenterol.* Nov 2007;102(11):2589-2600. PMID 17850410
- 9. Stebbing J, Farouk L, Panos G, et al. A meta-analysis of transient elastography for the detection of hepatic fibrosis. *J Clin Gastroenterol.* Mar 2010;44(3):214-219. PMID 19745758
- 10. Talwalkar JA, Kurtz DM, Schoenleber SJ, et al. Ultrasound-based transient elastography for the detection of hepatic fibrosis: systematic review and meta-analysis. *Clin Gastroenterol Hepatol.* Oct 2007;5(10):1214-1220. PMID 17916549
- 11. Tsochatzis EA, Gurusamy KS, Ntaoula S, et al. Elastography for the diagnosis of severity of fibrosis in chronic liver disease: a meta-analysis of diagnostic accuracy. *J Hepatol.* Apr 2011;54(4):650-659. PMID 21146892
- 12. Li Y, Huang YS, Wang ZZ, et al. Systematic review with meta-analysis: the diagnostic accuracy of transient elastography for the staging of liver fibrosis in patients with chronic hepatitis B. *Aliment Pharmacol Ther.* Feb 2016;43(4):458-469. PMID 26669632
- 13. Njei B, McCarty TR, Luk J, et al. Use of transient elastography in patients with HIV-HCV coinfection: A systematic review and meta-analysis. *J Gastroenterol Hepatol*. Oct 2016;31(10):1684-1693. PMID 26952020 14.Pavlov CS, Casazza G, Nikolova D, et al. Transient elastography for diagnosis of stages of hepatic fibrosis and cirrhosis in people with alcoholic liver disease. *Cochrane Database Syst Rev.* Jan 22 2015;1:CD010542. PMID 25612182
- 15.Xu X, Su Y, Song R, et al. Performance of transient elastography assessing fibrosis of single hepatitis B virus infection: a systematic review and meta-analysis of a diagnostic test. *Hepatol Int.* Oct 2015;9(4):558-566. PMID 26187292

---- CLICK THE ENVELOPE ICON BELOW TO SUBMIT COMMENTS

This medical policy is made available to you for informational purposes only. It is not a guarantee of payment or a substitute for your medical judgment in the treatment of your patients. Benefits and eligibility are determined by the member's subscriber agreement or member certificate and/or the employer agreement, and those documents will supersede the provisions of this medical policy. For information on member-specific benefits, call the provider call center. If you provide services to a member which are determined to not be medically necessary (or in some cases medically necessary services which are non-covered benefits), you may not charge the member for the services unless you have informed the member and they have agreed in writing in advance to continue with the treatment at their own expense. Please refer to your participation agreement(s) for the applicable provisions. This policy is current at the time of publication; however, medical practices, technology, and knowledge are constantly changing. BCBSRI reserves the right to review and revise this policy for any reason and at any time, with or without notice. Blue Cross & Blue Shield of Rhode Island is an independent licensee of the Blue Cross and Blue Shield Association.

